

SPECTRAGENIUS

- Release Notes -

CONTENT

-	<i>RELEASE NOTES –</i>	<i>1</i>
1.	<i>GENERAL INFORMATION</i>	<i>3</i>
2.	<i>SUPPORTED PLATFORMS</i>	<i>3</i>
3.	<i>KNOWN ISSUES</i>	<i>3</i>
4.	<i>WHAT TO DO WHEN YOU HAVE QUESTIONS.....</i>	<i>3</i>

1. General Information

This document contains information on the release of SpectraGenius 4.00 build 409.

SpectraGenius has been thoroughly tested. Some known issues are listed below.

2. Supported platforms

- SpectraGenius will work only on 64-bit Windows 8 and higher versions.
- A free USB-A port is required for the USB dongle.

3. Known issues

The following issues were identified. These are benign issues that can be worked around as described:

Function	Description of issue
Conversions & Corrections/Y Unit Conversions/to Log(1/R)	No undo available for conversion to Log(1/R). Workaround: Use the other function options to convert to the previous unit.
Data/Edit Title	Cannot undo/ Undo is only evident in the Edit Title screen, but does not apply to the Spectrum Area and File List. Workaround: After undoing, click “Apply” again.
Data/File Browser	Bruker file type filter is not applied in File Browser. Workaround: Use “Open with Windows Dialog”
Spectrum Display/left and right arrows	Please note that left and right are swapped (i.e., clicking the left arrow will move the spectrum to the right and vice versa) if the x axis unit is um or nm.

4. What to do when you have questions

If you have questions or experience difficulties with SpectraGenius, please do not hesitate to contact us for assistance:

support@stjapan.de or call us at +49 (0)2234 956372.